



BIDS AND AWARDS COMMITTEE

SUPPLEMENTAL BID BULLETIN

For ITB No. 2021-0207-ATQAL

PROCUREMENT OF SUPPLY DELIVERY, INSTALLATION, COMMISSIONING AND TESTING OF HIGH-PERFORMANCE LIQUID CHROMATOGRAPHS FOR ALABANG TESTING AND QUALITY ASSURANCE (ATQAL)

04 October 2021

In response to the issues/ concerns raised during the Pre-Bid Conference held on 30 September 2021 for the above mentioned ITB, hereunder are the amendments, clarifications and recommendations.

Refer to the comparative table given below for the previous (from) and revised (to) specifications:

A. Under Invitation to Bid for Submission of Price Quotation:

FROM	TO
Special Alloted General Funds (SAGF)	Special Account in the General Fund (SAGF)
Checklist of Documents (Legal Documents) Item (3) Statement of the bidder's Single Largest Completed Contract (SLCC) <u>similar to the contract to be bid</u> , except under conditions provided for in Sections 23.4.1.3 and 23.4.2.4 of the 2016 revised IRR of RA No. 9184, within the relevant period as provided in the Bidding Documents; and	Since the procurement allows partial bids, the procurement of which is divided into several lots, where each awarded lot gives rise to a separate contract, the SLCC requirement shall be based on the ABC for each lot (<u>GPPB NPM No. 149-2015</u>)
Payment of Bid Documents, if bid would be specific for a lot or line item	In case of lot procurement, the maximum fee for the bidding documents for each lot shall be based on its ABC, in accordance with the Guidelines issued by the GPPB



B. Under Section VII – Technical Specification:

DESCRIPTION	FROM	TO
ITEM NUMBERS: 1,2,3,4, and 5 (ALL ITEMS)		
A.3 COLUMN OVEN	Temperature control range $\geq -10^{\circ}\text{C}$ from room temperature or ambient temperature to 65°C	Temperature control range should include 20°C to 65°C
A.4 PHOTODIODE ARRAY DETECTOR	A.4 PHOTODIODE ARRAY DETECTOR	A.4 PHOTODIODE ARRAY DETECTOR or DIODE ARRAY DETECTOR
J. DELIVERY PERIOD	Not more than ninety (90) calendar days upon receipt of Notice to Proceed. To be delivered at FDA-ATQAL, Civic Drive, Filinvest Corporate City, Alabang, Muntinlupa City	Not more than one hundred twenty (120) calendar days upon receipt of Notice to Proceed. To be delivered at FDA-ATQAL, Civic Drive, Filinvest Corporate City, Alabang, Muntinlupa City
ITEM NUMBERS: 1,3, and 5 (HPLC-PDA; HPLC-PDA-ELD; UHPLC-PDA-ECD)		
C.4 STANDARDS (at least one (1) bottle each)	<p>NOTE:</p> <p><u>Traceability:</u> Standards must be traceable to NIST or equivalent international standard</p> <p><u>Expiration:</u> must be current lot for USP RS; for other standards, expiration, if possible, must not be less than one (1) year upon delivery.</p> <p><u>Certificate of analysis:</u> must be submitted upon delivery Presented below</p>	<p>NOTE:</p> <p><u>Traceability:</u> Standards must be traceable to NIST or equivalent international standard</p> <p><u>Expiration:</u> - must be current lot for USP RS; - <u>for other standards,</u> expiration, if possible, must not be less than one (1) year upon delivery. If one (1) year validity is not possible, please submit proof or supporting document/s upon delivery.</p> <p><u>Certificate of analysis:</u> must be submitted upon delivery Presented below</p>

ITEM NO. 1 SUPPLY, DELIVERY, INSTALLATION, COMMISSIONING AND TESTING OF HIGH PERFORMANCE LIQUID CHROMATOGRAPH WITH PHOTODIODE ARRAY DETECTOR (HPLC-PDA) COMPLETE WITH ALL NECESSARY PARTS AND ACCESSORIES FOR FDA-ALABANG TESTING QUALITY ASSURANCE LABORATORY (ATQAL)		
C.5 One (1) complete system of Filter Assembly	- 1 unit Vacuum pump (oil-free) - 1 pc. Borosilicate glass receiver (at least 4L) - 1 pc. Suction filter (funnel \geq 90mm)	- 1 unit Vacuum pump (oil-free) - 1 pc. Borosilicate glass receiver (at least 4L) - 1 pc. Suction filter (funnel minimum capacity: 1L)
C.7 OTHER LABORATORY SUPPLIES/ CONSUMABLES	a. At least 500 pcs. membrane filters (90mm diameter, 0.45 μ m pore size, hydrophilic PVDF)	a. At least 500 pcs. membrane filters (0.45 μ m pore size, hydrophilic PVDF). Note: diameter should be suitable to the suction filter funnel offered in C.5
I.5 The supplier must have at least ten (10) local installations of HPLC that can be viewed and demo to the end user, if necessary. Two (2) out of ten (10) local installations must be the same model as the one being offered.	The supplier must have at least ten (10) local installations of HPLC that can be viewed and demo to the end user, if necessary. Two (2) out of ten (10) local installations must be the same model as the one being offered.	The supplier must have at least ten (10) local installations of HPLC that can be viewed and demo to the end user, if necessary. Two (2) out of ten (10) local installations must be the same model and detector as the one being offered.
ITEM NO. 2 SUPPLY, DELIVERY, INSTALLATION, COMMISSIONING AND TESTING OF HIGH PERFORMANCE LIQUID CHROMATOGRAPH WITH PHOTODIODE ARRAY AND REFRACTIVE INDEX DETECTOR (HPLC-PDA-RID) COMPLETE WITH ALL NECESSARY PARTS AND ACCESSORIES FOR FDA-ALABANG TESTING QUALITY ASSURANCE LABORATORY (ATQAL)		
C.3 HPLC COLUMNS (1 pc. each)	a. L1 4.0mm x 250mm, 10 μ m b. L1 4.0mm x 300mm, 10 μ m c. L1 3.9mm x 300mm, 10 μ m	a. L1 4.0mm x 250mm, 10 μ m b. L1 4.0mm x 300mm, 5 μ m to 10 μ m c. L1 3.9mm x 300mm, 5 μ m to 10 μ m
C.5 One (1) complete system of Filter Assembly	C.5 One (1) complete system of Filter Assembly - 1 unit Vacuum pump (oil-free) - 1 pc. Borosilicate glass receiver (at least 4L) - 1 pc. Suction filter (funnel \geq 90mm)	C.4 One (1) complete system of Filter Assembly - 1 unit Vacuum pump (oil-free) - 1 pc. Borosilicate glass receiver (at least 4L) - 1 pc. Suction filter (funnel minimum capacity: 1L)


C.7 OTHER LABORATORY SUPPLIES/ CONSUMABLES	C.7 OTHER LABORATORY SUPPLIES/ CONSUMABLES a. At least 500 pcs. membrane filters (90mm diameter, 0.45µm pore size, hydrophilic PVDF)	C.6 OTHER LABORATORY SUPPLIES/ CONSUMABLES a. At least 500 pcs. membrane filters (0.45µm pore size, hydrophilic PVDF). Note: diameter should be suitable to the suction filter funnel offered in C.4 (One (1) complete system of Filter Assembly)
CORRECTION ON NUMBERING	- C.6 - C.7 - C.8 - C.9 - C.10 - C.11	- C.5 - C.6 - C.7 - C.8 - C.9 - C.10
ITEM NO. 3 SUPPLY, DELIVERY, INSTALLATION, COMMISSIONING AND TESTING OF HIGH PERFORMANCE LIQUID CHROMATOGRAPH WITH PHOTODIODE ARRAY AND FLUORESCENCE DETECTOR (HPLC-PDA-FLD) COMPLETE WITH ALL NECESSARY PARTS AND ACCESSORIES FOR FDA-ALABANG TESTING QUALITY ASSURANCE LABORATORY (ATQAL)		
A.5 MULTIWAVELENGTH FLUORESCENCE DETECTOR	A.5 MULTIWAVELENGTH FLUORESCENCE DETECTOR	A.5 MULTIWAVELENGTH FLUORESCENCE DETECTOR OR FLUORESCENCE DETECTOR
C.5 One (1) complete system of Filter Assembly	- 1 unit Vacuum pump (oil-free) - 1 pc. Borosilicate glass receiver (at least 4L) - 1 pc. Suction filter (funnel ≥ 90mm)	- 1 unit Vacuum pump (oil-free) - 1 pc. Borosilicate glass receiver (at least 4L) - 1 pc. Suction filter (funnel minimum capacity: 1L)
C.7 OTHER LABORATORY SUPPLIES/ CONSUMABLES	a. At least 500 pcs. membrane filters (90mm diameter, 0.45µm pore size, hydrophilic PVDF)	a. At least 500 pcs. membrane filters (0.45µm pore size, hydrophilic PVDF). Note: diameter should be suitable to the suction filter funnel offered in C.5
	b. At least 100 pcs. Magnetic Stir Bar without pivot ring (25mm x 5mm)	At least 100 pcs. Magnetic Stir Bar without pivot ring (25mm x 5mm or 25mm x 8mm)

ITEM NO. 4 SUPPLY, DELIVERY, INSTALLATION, COMMISSIONING AND TESTING OF HIGH PERFORMANCE LIQUID CHROMATOGRAPH WITH PHOTODIODE ARRAY AND UV-VIS DETECTOR (HPLC-PDA-UV-VIS) COMPLETE WITH ALL NECESSARY PARTS AND ACCESSORIES FOR FDA-ALABANG TESTING QUALITY ASSURANCE LABORATORY (ATQAL)

I.5	The supplier must have at least ten (10) local installations of HPLC that can be viewed and demo to the end user, if necessary. Two (2) out of ten (10) local installations must be the same model as the one being offered.	The supplier must have at least ten (10) local installations of HPLC that can be viewed and demo to the end user, if necessary. Two (2) out of ten (10) local installations must be the same model with either UV-VIS or PDA detector , as the one being offered.
CORRECTION ON NUMBERING	- C.7	- C.4

ITEM NO. 5 SUPPLY, DELIVERY, INSTALLATION, COMMISSIONING AND TESTING OF ULTRA HIGH PERFORMANCE LIQUID CHROMATOGRAPH WITH PHOTODIODE ARRAY AND ELECTROCHEMICAL DETECTOR (UHPLC-PDA-ECD) COMPLETE WITH ALL NECESSARY PARTS AND ACCESSORIES FOR FDA-ALABANG TESTING QUALITY ASSURANCE LABORATORY (ATQAL)

I.5	<p>The supplier must have at least ten (10) local installations of HPLC that can be viewed and demo to the end user, if necessary. Two (2) out of ten (10) local installations must be the same model as the one being offered.</p> <p>Note: Upon submission of bid, the bidder must submit a list of companies with installed unit. Contact person and contact details (telephone number/mobile number and e-mail address) for each listed company must also be submitted.</p>	<p>The supplier must have at least ten (10) local or within Asia-Pacific installations of UHPLC-PDA-ECD or equivalent UHPLC with either PDA and/or ECD detectors that can be viewed and demo to the end user, if necessary.</p> <p>Note: Upon submission of bid, the bidder must submit a list of companies with installed unit. Contact person and contact details (telephone number/mobile number and e-mail address) for each listed company must also be submitted.</p>
-----	--	---


ENGR. ANA TRINIDAD F. RIVERA, MSc
 Chairperson, FDA-Bids and Awards Committee