

Description	Packaging	QTY	UNIT COST
Anaerobic Jar, 3.5L Hold up to 15 petridish	Piece	1	22,000.00
Autoclavable Specimen Bags Disposable, polypropylene, approximately 1.5 mm thick, Dimensions: 55 cm x 75 cm (W x H) or bigger	200 pieces/box	1	13,000.00
Autoclavable Polypropylene Test Tube Closure, 20mm	250 pieces/pack	1	5,400.00
Biohazard bag stand/holder with cover, holds 21.59 cm x 27.94 cm bag, autoclavable	piece	25	2,000.00
Biohazard Disposable Bags, 35.56 cm x 48.26 cm, with Biohazard print	100 pieces/pack	2	6,650.00
Biohazard Disposable Bags, small, 21.59 cm x 27.94 cm, Polypropylene (PP) with Biohazard print	100 pieces/pack	3	2,700.00
Cooling Block for 0.5mL Tubes, Strips or 96-Well Plates	piece	1	39,000.00
Cooling Block for 1.5-2mL tubes, 24-well	piece	1	27,000.00
Cryopreservation beads, 20 beads in 8 ml vial	64 pieces/box	1	3,000.00
Cryovials 1.8 mL, nuclease-free Slim type (RNA and DNA free)	500pieces/pack	1	20,000.00
Low-profile 0.2 ml 8-strip tube without caps	120 strip/pack	2	5,000.00
Microcentrifuge tube opener (Cap Opener for 0.5, 0.65, 1.5, 1.7 and 2.0mL snap cap microcentrifuge tubes), Autoclavable, Polypropylene	5 pieces/ Case	2	5,000.00
PCR Box/Rack, Holds 96-well PCR plates, tube strips, and individual 0.2mL tubes	5pieces/pack	2	5,000.00
Petri Dish Can, Holds 10 dishes 100 mm x 15 mm, staineless steel, with removable rack	piece	2	5,000.00
Pipette Tips, 0.1 mL, Autoclavable, compatible with BRAND Transferpette Pipettors Note: sample to be submitted for end-users verification	pack of 200	1	3,500.00
Pipette Tips, 0.5-10µL, disposable, compatible with BRAND Transferpette Pipettors Note: sample tip should be provided for end- user/TWG evaluation	200 pieces/pack	6	1,000.00

Pipette Tips, 0.5-5mL,disposable, compatible withTrefflab Transferpette pipettors Note: sample tip should be provided for end- user/TWG evaluation	200 pieces/pack	3	3,500.00
Pipette Tips, 1.0 mL, Autoclavable, compatible with BRAND Transferpette or Volac Pipettors Note: sample to be submitted for end-users verification	200 pieces/ pack	1	2,000.00
Pipette Tips, 10 mL, Autoclavable, compatible with BRAND Transferpette pipettor Note: sample tip should be provided for end- user/TWG evaluation	200 pieces/ pack	2	3,000.00
Pipette Tips, 100-1000µL, disposable, compatible with BRAND Transferpette pipettor Note: sample tip should be provided for end- user/TWG evaluation	200 pieces/pack	8	2,000.00
Pipette Tips, filtered, Sterile 20-200 µL, disposable, compatible with BRAND Transferpette pipettor Note: sample tip should be provided for end- user/TWG evaluation	96 pieces /rack, 10 racks/unit	2	7,500.00
Pipette Tips, filtered, Sterile 30-300 µL, disposable, compatible with BRAND Transferpette pipettor Note: sample tip should be provided for end- user/TWG evaluation	96 pieces /rack, 10 racks/unit	2	8,500.00
Plate holder, 96 wells, holds twelve 1 x 8 strip tubes	20 pieces/case	1	8,000.00
Racked Sterile Filter Pipette Tips, 100-1000 uL	960/case	2	4,900.00
Reagent reservior (polypropylene, autoclavable, 60 mL capacity)	5 pieces/pack	10	1,000.00

<p>Reversible Racks for Microtubes 80/60 well rack capacities : 80-place for 1.5mL and 2mL tubes; 60-place for 0.5mL and 0.2 mL tubes Dimensions (L x W x H): 80/60 Well: 228 x 67 x 27mm (9" x 2.6" x 1") Wells are easily identified with molded numbers Made with heavy duty polypropylene Ideal for freezing and autoclaving</p>	5 pieces/ pack	3	5,000.00
<p>Reversible racks for microtubes 96 well rack capacities: 96 places for 0.5 mL microtubes on one side or 1.5, 1.7, 2.0 mL microtubes on reverse side Wells are easily identified with molded numbers Dimensions (L x W x H): 221 x 131 x 46mm (8.7" x 5.2" x 1.8") Made of heavy duty polypropylene Lid/cover included Ideal for freezing and autoclaving</p>	5 pieces/ pack	2	8,000.00
<p>Snaplock microtubes, 1.5ml, nuclease-free</p>	500 pieces/pack	2	6,000.00
<p>Sterility test canister, dual-needle metal spike, sterile, vented, ideal for sampling for antibiotic products 0.45µm membrane,</p>	10 pieces/box	2	48,000.00