

Description	Packaging	QTY	UNIT COST
Anaerobic Jar, 3.5L Hold up to 15 petridish	Piece	1	22,000.00
Autoclavable Specimen Bags Disposable, polypropylene, approximately 1.5 mm thick, Dimensions: 55 cm x 75 cm (W x H) or bigger	200 pieces/box	1	13,000.00
Pipette rack with 5mL pipette tips, autocleavable, compatible with BRAND Socorex Acura 835 Note: sample tip should be provided for end-user/TWG	at least 24 tips/rack	3	4,000.00
Pipette Tips, 0-5mL, Autoclavable, compatible with BRAND Socorex Acura 835 Note: sample tip should be provided for end-user/TWG evaluation	200 pieces/ pack	3	5,000.00
Pipette Tips, 0.5-10µL, disposable, compatible with BRAND Transferpette Pipettors Note: sample tip should be provided for end-user/TWG evaluation	200 pieces/pack	1	1,000.00
Pipette tips, 0.5-200µL, compatible with Micropipette Smart Accumax (20-200uL) Note: sample tip should be provided for end- user/TWG	1000 pieces/pack	1	760.00
Pipette Tips, 100-1000µL, disposable, compatible with BRAND Transferpette pipettor Note: sample tip should be provided for end-user/TWG evaluation	200 pieces/pack	1	2,000.00
Pipette Tips, blue, 100uL - 1000uL, Autoclavable, compatible with BRAND Transferpette Pipettor and/or Nichiryo Pipettor Note: sample tip should be provided for end-user/TWG evaluation	pack of 1000	5	1,650.00
Sample bag with full-page filter, 190 mm X 300mm, maximum blending volume up to 400 mL, sterile, designed for stomacher blenders	500 pieces/box	1	21,000.00
Soldering Machine Electric Soldering Iron Kit with soldering tip holder, handheld, 220V, 100W, 50Hz/60Hz, adjustable temperature with LCD, with different soldering iron tips and consummables Temperature range: 180-500 degree celcius	piece	1	5,000.00
Sterility test canister, dual-needle metal spike, sterile, vented, ideal for sampling LVP style containers. 0.45um membrane. (10 pcs/box)	10 pieces/box	4	48,000.00