

FDA PRAISE Loyalty Incentive – Watch tag



TECHNICAL SPECIFICATION

Republic of the Philippines
Food and Drug Administration

Lot No. 1 Item No. 1	Watch -Female for 10 years in government service	Qty. / Unit	12
---------------------------------------	---	--------------------	-----------

Name of Manufacturer:	Country of Origin:
------------------------------	---------------------------

Brand:	Model: (if applicable)
---------------	-------------------------------

Total ABC for Lot No. 1:**PURCHASER'S SPECIFICATION****STATEMENT OF COMPLIANCE****A. Detailed Technical Specifications:**

1. Case size: at least 32mm
2. Case material: Stainless steel
3. Glass Material: Mineral Glass
4. Strap Materials: Leather; black
5. Dial, Case & Strap Color: Must be complementary (white/silver/black)
6. Movement: Quartz
7. Dial Hands: Basic 3 hands
8. Water Resistant
9. Scratch Resistant
10. Watch Closure: Buckle
11. Branded
12. Warranty: 1 year

B. Additional Documentary Requirements:

1. The bidder shall submit any of the following whichever is applicable:
 - a. If the bidder is a manufacturer, certify that the bidder manufactures the products/items; or
 - b. If the bidder is an agent of the exclusive distributor or dealer, the following must be provided:
 - i. Certificate or Distributor/Dearlership Agreement by the Manufacturer with the distributor or dealer; and
 - ii. Contract between the distributor/dealer and the bidder
2. Warranty certificate

C. Instruction before mass production: Sample of watch for each year subject for approval of the End-user before mass production.**D. Upon delivery the following shall be complied:**

1. Each watch must have an attached tag -- silver plated tag with embossed/engraved FDA logo and no. of years in service
- 2.. Packaging Instructions:
 - a. Primary Packaging: Each piece is individually place in a color light gray PU leather box with pillow
 - b. Secondary Packaging: individually packed, watched must be placed in a bigger box with the following label:
 - i. Watch for: specify gender and no. of years in service
 - ii. Total no. of watch per box

TECHNICAL SPECIFICATION			
Republic of the Philippines Food and Drug Administration			
Lot No. 1	Watch -Female for 15 years in government service	Qty. / Unit	6
Item No.			
Name of Manufacturer:		Country of Origin:	
Brand:		Model: (if applicable)	
Total ABC for Lot No. 1:			
PURCHASER'S SPECIFICATION		STATEMENT OF COMPLIANCE	
A. Detailed Technical Specifications: 1. Case size: at least 35mm 2. Case material: Stainless steel 3. Glass Material: Mineral Glass 4. Strap Materials: Leather; black 5. Dial, Case & Strap Color: Must be complementary (white/silver/black) 6. Movement: Quartz 7. Dial Hands: Basic 3 hands 8. Water Resistant 9. Scratch Resistant 10. Watch Closure: Buckle 11. Branded 12. Warranty: 1 year			
B. Additional Documentary Requirements: 1. The bidder shall submit any of the following whichever is applicable: a. If the bidder is a manufacturer, certify that the bidder manufactures the products/items; or b. If the bidder is an agent of the exclusive distributor or dealer, the following must be provided: i. Certificate or Distributor/Dearlership Agreement by the Manufacturer with the distributor or dealer; and ii. Contract between the distributor/dealer and the bidder 2. Warranty certificate			
C. Instruction before mass production: Sample of watch for each year subject for approval of the End-user before mass production.			
D. Upon delivery the following shall be complied: 1. Each watch must have an attached tag -- silver plated tag with embossed/engraved FDA logo and no. of years in service 2.. Packaging Instructions: a. Primary Packaging: Each piece is individually place in a color light gray PU leather box with pillow b. Secondary Packaging: individually packed, watched must be placed in a bigger box with the following label: i. Watch for: specify gender and no. of years in service ii. Total no. of watch per box			

TECHNICAL SPECIFICATION

Republic of the Philippines
Food and Drug Administration

Lot No. 1	Watch -Female for 20 years in government service	Qty. / Unit	12
Item No.			

Name of Manufacturer:	Country of Origin:
------------------------------	---------------------------

Brand:	Model: (if applicable)
---------------	-------------------------------

Total ABC for Lot No. 1:

PURCHASER'S SPECIFICATION	STATEMENT OF COMPLIANCE
---------------------------	-------------------------

A. Detailed Technical Specifications:

1. Case size: at least 30mm
2. Case material: Stainless steel
3. Glass Material: Mineral Glass
4. Strap Materials: Silver stainless steel; alloy band
5. Dial, Case & Strap Color: Must be complementary (white/silver/black)
6. Movement: Japanese Quartz
7. Water Resistant
8. Scratch Resistant
9. Watch Closure: Clasp
10. Branded
11. Warranty: 1 year

B. Additional Documentary Requirements:

1. The bidder shall submit any of the following whichever is applicable:
 - a. If the bidder is a manufacturer, certify that the bidder manufactures the products/items; or
 - b. If the bidder is an agent of the exclusive distributor or dealer, the following must be provided:
 - i. Certificate or Distributor/Dearlership Agreement by the Manufacturer with the distributor or dealer; and
 - ii. Contract between the distributor/dealer and the bidder
2. Warranty certificate

C. Instruction before mass production: Sample of watch for each year subject for approval of the End-user before mass production.

D. Upon delivery the following shall be complied:

1. Each watch must have an attached tag -- silver plated tag with embossed/engraved FDA logo and no. of years in service
2. Packaging Instructions:
 - a. Primary Packaging: Each piece is individually placed in a color light gray PU leather box with pillow
 - b. Secondary Packaging: individually packed, watched must be placed in a bigger box with the following label:
 - i. Watch for: specify gender and no. of years in service
 - ii. Total no. of watch per box

TECHNICAL SPECIFICATION			
Republic of the Philippines Food and Drug Administration			
Lot No. 1 Item No.	Watch -Female for 25 years in government service	Qty. / Unit	16
Name of Manufacturer:		Country of Origin:	
Brand:		Model: (if applicable)	
Total ABC for Lot No. 1:			
PURCHASER'S SPECIFICATION		STATEMENT OF COMPLIANCE	
A. Detailed Technical Specifications: 1. Case size: at least 36mm 2. Case thickness: 6-8mm 3. Case material: Stainless steel 4. Glass Material: Mineral Glass 5. Strap Materials: Silver stainless steel; width: at least 14mm 6. Dial, Case & Strap Color: Must be complementary (silver) 7. Movement: Japanese Quartz 8. Water Resistant 9. Scratch Resistant 10. Watch Closure: Double Press Deployant 11. Branded 12. Warranty: 1 year			
B. Additional Documentary Requirements: 1. The bidder shall submit any of the following whichever is applicable: a. If the bidder is a manufacturer, certify that the bidder manufactures the products/items; or b. If the bidder is an agent of the exclusive distributor or dealer, the following must b provided: i. Certificate or Distributor/Dearlanship Agreement by the Manufacturer with the distributor or dealer; and ii. Contract between the distributor/dealer and the bidder 2. Warranty certificate			
C. Instruction before mass production: Sample of watch for each year subject for approval of the End-user before mass production.			
D. Upon delivery the following shall be complied: 1. Each watch must have an attached tag -- silver plated tag with embossed/engraved FDA logo and no. of years in service 2.. Packaging Instructions: a. Primary Packaging: Each piece is individually place in a color light gray PU leather box with pillow b. Secondary Packaging: individually packed, watched must be placed in a bigger box with the following label: i. Watch for: specify gender and no. of years in service ii. Total no. of watch per box			

TECHNICAL SPECIFICATION			
Republic of the Philippines Food and Drug Administration			
Lot No. 1	Watch -Female for 30 years in government service	Qty. / Unit	21
Item No.			
Name of Manufacturer:		Country of Origin:	
Brand:		Model: (if applicable)	
Total ABC for Lot No. 1:			
PURCHASER'S SPECIFICATION		STATEMENT OF COMPLIANCE	
A. Detailed Technical Specifications: 1. Case size: at least 35mm 2. Case thickness: 8-10mm 3. Case material: Stainless steel 4. Glass Material: Mineral Glass with sapphire 5. Strap Materials: Silver stainless steel; width: at least 16mm 6. Dial, Case & Strap Color: Must be complementary (Gold) 7. Movement: Quartz Multifunction 8. Water Resistant 9. Scratch Resistant 10. Watch Closure: Double Press Deployant 11. Branded 12. Warranty: 1 year			
B. Additional Documentary Requirements: 1. The bidder shall submit any of the following whichever is applicable: a. If the bidder is a manufacturer, certify that the bidder manufactures the products/items; or b. If the bidder is an agent of the exclusive distributor or dealer, the following must b provided: i. Certificate or Distributor/Dearlership Agreement by the Manufacturer with the distributor or dealer; and ii. Contract between the distributor/dealer and the bidder 2. Warranty certificate			
C. Instruction before mass production: Sample of watch for each year subject for approval of the End-user before mass production.			
D. Upon delivery the following shall be complied: 1. Each watch must have an attached tag -- silver plated tag with embossed/engraved FDA logo and no. of years in service 2.. Packaging Instructions: a. Primary Packaging: Each piece is individually place in a color light gray PU leather box with pillow b. Secondary Packaging: individually packed, watched must be placed in a bigger box with the following label: i. Watch for: specify gender and no. of years in service ii. Total no. of watch per box			

TECHNICAL SPECIFICATION			
Republic of the Philippines Food and Drug Administration			
Lot No. 1	Watch -Female for 35 years in government service	Qty. / Unit	7
Item No.			
Name of Manufacturer:		Country of Origin:	
Brand:		Model: (if applicable)	
Total ABC for Lot No. 1:			
PURCHASER'S SPECIFICATION		STATEMENT OF COMPLIANCE	
A. Detailed Technical Specifications: 1. Case size: at least 36mm 2. Case thickness: 8-10mm 3. Case material: Stainless steel 4. Glass Material: Mineral Glass with sapphire 5. Strap Materials: stainless steel; width: at least 16mm 6. Dial, Case & Strap Color: Must be complementary (Gold) 7. Movement: Quartz Multifunction 8. Water Resistant 9. Scratch Resistant 10. Watch Closure: Fold over clasp with push button 11. Branded 12. Warranty: 1 year			
B. Additional Documentary Requirements: 1. The bidder shall submit any of the following whichever is applicable: a. If the bidder is a manufacturer, certify that the bidder manufactures the products/items; or b. If the bidder is an agent of the exclusive distributor or dealer, the following must b provided: i. Certificate or Distributor/Dearlership Agreement by the Manufacturer with the distributor or dealer; and ii. Contract between the distributor/dealer and the bidder 2. Warranty certificate			
C. Instruction before mass production: Sample of watch for each year subject for approval of the End-user before mass production.			
D. Upon delivery the following shall be complied: 1. Each watch must have an attached tag -- silver plated tag with embossed/engraved FDA logo and no. of years in service 2.. Packaging Instructions: a. Primary Packaging: Each piece is individually place in a color light gray PU leather box with pillow b. Secondary Packaging: individually packed, watched must be placed in a bigger box with the following label: i. Watch for: specify gender and no. of years in service ii. Total no. of watch per box			

TECHNICAL SPECIFICATION			
Republic of the Philippines Food and Drug Administration			
Lot No. 1 Item No.	Watch -Female for 40 years in government service	Qty. / Unit	8
Name of Manufacturer:		Country of Origin:	
Brand:		Model: (if applicable)	
Total ABC for Lot No. 1:			
PURCHASER'S SPECIFICATION		STATEMENT OF COMPLIANCE	
A. Detailed Technical Specifications: 1. Case size: at least 36mm 2. Case thickness: 8-10mm 3. Case material: Stainless steel 4. Glass Material: Mineral Glass with sapphire 5. Strap Materials: stainless steel; width: at least 16mm 6. Dial, Case & Strap Color: Must be complementary (Gold) 7. Movement: Automatic 8. Water Resistant 9. Scratch Resistant 10. Watch Closure: Safety mesh buckle 11. Branded 12. Warranty: 1 year			
B. Additional Documentary Requirements: 1. The bidder shall submit any of the following whichever is applicable: a. If the bidder is a manufacturer, certify that the bidder manufactures the products/items; or b. If the bidder is an agent of the exclusive distributor or dealer, the following must b provided: i. Certificate or Distributor/Dearlership Agreement by the Manufacturer with the distributor or dealer; and ii. Contract between the distributor/dealer and the bidder 2. Warranty certificate			
C. Instruction before mass production: Sample of watch for each year subject for approval of the End-user before mass production.			
D. Upon delivery the following shall be complied: 1. Each watch must have an attached tag -- silver plated tag with embossed/engraved FDA logo and no. of years in service 2.. Packaging Instructions: a. Primary Packaging: Each piece is individually place in a color light gray PU leather box with pillow b. Secondary Packaging: individually packed, watched must be placed in a bigger box with the following label: i. Watch for: specify gender and no. of years in service ii. Total no. of watch per box			

TECHNICAL SPECIFICATION			
Republic of the Philippines Food and Drug Administration			
Lot No. 1	Watch -Male for 10 years in government service	Qty. / Unit	6
Item No.			
Name of Manufacturer:		Country of Origin:	
Brand:		Model: (if applicable)	
Total ABC for Lot No. 1:			
PURCHASER'S SPECIFICATION		STATEMENT OF COMPLIANCE	
A. Detailed Technical Specifications: 1. Case size: at least 36mm 2. Case material: Stainless steel 3. Glass Material: Mineral Glass 4. Strap Materials: Leather; black 5. Dial, Case & Strap Color: Must be complementary (white/silver/black) 6. Movement: Quartz 7. Water Resistant 8. Watch Closure: Pin buckle fastening 9. Branded 10. Warranty: 1 year			
B. Additional Documentary Requirements: 1. The bidder shall submit any of the following whichever is applicable: a. If the bidder is a manufacturer, certify that the bidder manufactures the products/items; or b. If the bidder is an agent of the exclusive distributor or dealer, the following must b provided: i. Certificate or Distributor/Dearlanship Agreement by the Manufacturer with the distributor or dealer; and ii. Contract between the distributor/dealer and the bidder 2. Warranty certificate C. Instruction before mass production: Sample of watch for each year subject for approval of the End-user before mass production.			
D. Upon delivery the following shall be complied: 1. Each watch must have an attached tag -- silver plated tag with embossed/engraved FDA logo and no. of years in service 2.. Packaging Instructions: a. Primary Packaging: Each piece is individually place in a color light gray PU leather box with pillow b. Secondary Packaging: individually packed, watched must be placed in a bigger box with the following label: i. Watch for: specify gender and no. of years in service ii. Total no. of watch per box			

TECHNICAL SPECIFICATION			
Republic of the Philippines Food and Drug Administration			
Lot No. 1	Watch -Male for 20 years in government service	Qty. / Unit	7
Item No.			
Name of Manufacturer:		Country of Origin:	
Brand:		Model: (if applicable)	
Total ABC for Lot No. 1:			
PURCHASER'S SPECIFICATION		STATEMENT OF COMPLIANCE	
A. Detailed Technical Specifications: 1. Case width: at least 40mm 2. Case thickness: at least 10mm 3. Case material: Stainless steel 4. Glass Material: Mineral sapphire glass 5. Strap Materials: Stainless steel; width at least 22mm 6. Dial, Case & Strap Color: Must be complementary (silver/black/white) 7. Movement: Quartz 8. Water Resistant 9. Watch Closure: Single Pusher Fold over Clasp 10. Branded 11. Warranty: 1 year			
B. Additional Documentary Requirements: 1. The bidder shall submit any of the following whichever is applicable: a. If the bidder is a manufacturer, certify that the bidder manufactures the products/items; or b. If the bidder is an agent of the exclusive distributor or dealer, the following must be provided: i. Certificate or Distributor/Dearlership Agreement by the Manufacturer with the distributor or dealer; and ii. Contract between the distributor/dealer and the bidder 2. Warranty certificate			
C. Instruction before mass production: Sample of watch for each year subject for approval of the End-user before mass production.			
D. Upon delivery the following shall be complied: 1. Each watch must have an attached tag -- silver plated tag with embossed/engraved FDA logo and no. of years in service 2.. Packaging Instructions: a. Primary Packaging: Each piece is individually place in a color light gray PU leather box with pillow b. Secondary Packaging: individually packed, watched must be placed in a bigger box with the following label: i. Watch for: specify gender and no. of years in service ii. Total no. of watch per box			

TECHNICAL SPECIFICATION			
Republic of the Philippines Food and Drug Administration			
Lot No. 1 Item No.	Watch -Male for 25 years in government service	Qty. / Unit	3
Name of Manufacturer:		Country of Origin:	
Brand:		Model: (if applicable)	
Total ABC for Lot No. 1:			
PURCHASER'S SPECIFICATION		STATEMENT OF COMPLIANCE	
<p>A. Detailed Technical Specifications:</p> <ol style="list-style-type: none"> 1. Case size: at least 42mm 2. Case thickness: at least 11mm 3. Case material: Stainless steel 4. Glass Material: Mineral sapphire glass 5. Strap Materials: Silver stainless steel; width at least 22mm 6. Dial, Case & Strap Color: Must be complementary (white/silver/black) 7. Movement: Quartz 8. Water Resistant 9. Watch Closure: 2 Pusher Fold over Clasp 10. Branded 11. Warranty: 1 year 			
<p>B. Additional Documentary Requirements:</p> <ol style="list-style-type: none"> 1. The bidder shall submit any of the following whichever is applicable: <ol style="list-style-type: none"> a. If the bidder is a manufacturer, certify that the bidder manufactures the products/items; or b. If the bidder is an agent of the exclusive distributor or dealer, the following must b provided: <ol style="list-style-type: none"> i. Certificate or Distributor/Dearlership Agreement by the Manufacturer with the distributor or dealer; and ii. Contract between the distributor/dealer and the bidder 			
<p>2. Warranty certificate</p>			
<p>C. Instruction before mass production: Sample of watch for each year subject for approval of the End-user before mass production.</p>			
<p>D. Upon delivery the following shall be complied:</p> <ol style="list-style-type: none"> 1. Each watch must have an attached tag -- silver plated tag with embossed/engraved FDA logo and no. of years in service 2.. Packaging Instructions: <ol style="list-style-type: none"> a. Primary Packaging: Each piece is individually place in a color light gray PU leather box with pillow b. Secondary Packaging: individually packed, watched must be placed in a bigger box with the following label: <ol style="list-style-type: none"> i. Watch for: specify gender and no. of years in service ii. Total no. of watch per box 			

TECHNICAL SPECIFICATION			
Republic of the Philippines Food and Drug Administration			
Lot No. 1 Item No.	Watch -Male for 30 years in government service	Qty. / Unit	8
Name of Manufacturer:		Country of Origin:	
Brand:		Model: (if applicable)	
Total ABC for Lot No. 1:			
PURCHASER'S SPECIFICATION		STATEMENT OF COMPLIANCE	
<p>A. Detailed Technical Specifications:</p> <ol style="list-style-type: none"> 1. Case size: at least 44mm 2. Case thickness: at least 8mm 3. Case material: Stainless steel 4. Glass Material: Mineral sapphire glass 5. Strap Materials: Silver stainless steel; width at least 22mm 6. Dial, Case & Strap Color: Must be complementary (Gold) 7. Movement: Quartz 8. Water Resistant 9. Watch Closure: 2 Pusher Fold over Clasp 10. Branded 11. Warranty: 1 year 			
<p>B. Additional Documentary Requirements:</p> <ol style="list-style-type: none"> 1. The bidder shall submit any of the following whichever is applicable: <ol style="list-style-type: none"> a. If the bidder is a manufacturer, certify that the bidder manufactures the products/items; or b. If the bidder is an agent of the exclusive distributor or dealer, the following must b provided: <ol style="list-style-type: none"> i. Certificate or Distributor/Dearlership Agreement by the Manufacturer with the distributor or dealer; and ii. Contract between the distributor/dealer and the bidder 2. Warranty certificate 			
<p>C. Instruction before mass production: Sample of watch for each year subject for approval of the End-user before mass production.</p>			
<p>D. Upon delivery the following shall be complied:</p> <ol style="list-style-type: none"> 1. Each watch must have an attached tag -- silver plated tag with embossed/engraved FDA logo and no. of years in service 2.. Packaging Instructions: <ol style="list-style-type: none"> a. Primary Packaging: Each piece is individually place in a color light gray PU leather box with pillow b. Secondary Packaging: individually packed, watched must be placed in a bigger box with the following label: <ol style="list-style-type: none"> i. Watch for: specify gender and no. of years in service ii. Total no. of watch per box 			

TECHNICAL SPECIFICATION

**Republic of the Philippines
Food and Drug Administration**

Lot No. 1	Watch -Male for 40 years in government service	Qty. / Unit	4
Item No.			

Name of Manufacturer:	Country of Origin:
------------------------------	---------------------------

Brand:	Model: (if applicable)
---------------	-------------------------------

Total ABC for Lot No. 1:

PURCHASER'S SPECIFICATION

STATEMENT OF COMPLIANCE

A. Detailed Technical Specifications:

1. Case size: 44-45mm
2. Case thickness: at least 8mm
3. Case material: Stainless steel
4. Glass Material: Mineral sapphire glass
5. Strap Materials: Silver stainless steel; width at least 22mm
6. Dial, Case & Strap Color: Must be complementary (Gold)
7. Movement: Automatic
8. Water Resistant
9. Watch Closure: 2 Pusher Fold over Clasp
10. CalendarType: Day-date
11. Features: Chronograph
12. Branded
13. Warranty: 1 year

B. Additional Documentary Requirements:

1. The bidder shall submit any of the following whichever is applicable:
 - a. If the bidder is a manufacturer, certify that the bidder manufactures the products/items; or
 - b. If the bidder is an agent of the exclusive distributor or dealer, the following must b provided:
 - i. Certificate or Distributor/Dearlership Agreement by the Manufacturer with the distributor or dealer; and
 - ii. Contract between the distributor/dealer and the bidder

2. Warranty certificate

C. Instruction before mass production: Sample of watch for each year subject for approval of the End-user before mass production.

D. Upon delivery the following shall be complied:

1. Each watch must have an attached tag -- silver plated tag with embossed/engraved FDA logo and no. of years in service
- 2.. Packaging Instructions:
 - a. Primary Packaging: Each piece is individually place in a color light gray PU leather box with pillow
 - b. Secondary Packaging: individually packed, watched must be placed in a bigger box with the following label:
 - i. Watch for: specify gender and no. of years in service
 - ii. Total no. of watch per box