



**TECHNICAL SPECIFICATION**

Republic of the Philippines  
Food and Drug Administration

Lot No. 1	Watch -Male for 40 years in government service	Qty. / Unit	4
Item No.			

Name of Manufacturer:	Country of Origin:
-----------------------	--------------------

Brand:	Model: (if applicable)
--------	------------------------

**Total ABC for Lot No. 1:****PURCHASER'S SPECIFICATION****STATEMENT OF COMPLIANCE****A. Detailed Technical Specifications:**

1. Case size: 44-45mm
2. Case thickness: at least 8mm
3. Case material: Stainless steel
4. Glass Material: Mineral sapphire glass
5. Strap Materials: Silver stainless steel; width at least 22mm
6. Dial, Case & Strap Color: Must be complementary (Gold)
7. Movement: Automatic
8. Water Resistant
9. Watch Closure: 2 Pusher Fold over Clasp
10. CalendarType: Day-date
11. Features: Chronograph
12. Branded
13. Warranty: 1 year

**B. Upon delivery the following shall be complied:**

1. Each watch must have an attached tag -- silver plated tag with embossed/engraved FDA logo and no. of years in service
- 2.. Packaging Instructions:
  - a. Primary Packaging: Each piece is individually place in a color light gray PU leather box with pillow
  - b. Secondary Packaging: individually packed, watched must be placed in a bigger box with the following label:
    - i. Watch for: specify gender and no. of years in service
    - ii. Total no. of watch per box

TECHNICAL SPECIFICATION			
Republic of the Philippines Food and Drug Administration			
Lot No. 1 Item No.	Watch -Female for 40 years in government service	Qty. / Unit	8
Name of Manufacturer:		Country of Origin:	
Brand:		Model: (if applicable)	
Total ABC for Lot No. 1:			
PURCHASER'S SPECIFICATION		STATEMENT OF COMPLIANCE	
<p><b>A. Detailed Technical Specifications:</b></p> <ol style="list-style-type: none"> <li>1. Case size: at least 36mm</li> <li>2. Case thickness: 8-10mm</li> <li>3. Case material: Stainless steel</li> <li>4. Glass Material: Mineral Glass with sapphire</li> <li>5. Strap Materials: stainless steel; width: at least 16mm</li> <li>6. Dial, Case &amp; Strap Color: Must be complementary (Gold)</li> <li>7. Movement: Automatic</li> <li>8. Water Resistant</li> <li>9. Scratch Resistant</li> <li>10. Watch Closure: Safety mesh buckle</li> <li>11. Branded</li> <li>12. Warranty: 1 year</li> </ol>			
<p><b>B. Upon delivery the following shall be complied:</b></p> <ol style="list-style-type: none"> <li>1. Each watch must have an attached tag -- silver plated tag with embossed/engraved FDA logo and no. of years in service</li> <li>2.. Packaging Instructions:               <ol style="list-style-type: none"> <li>a. Primary Packaging: Each piece is individually place in a color light gray PU leather box with pillow</li> <li>b. Secondary Packaging: individually packed, watched must be placed in a bigger box with the following label:                   <ol style="list-style-type: none"> <li>i. Watch for: specify gender and no. of years in service</li> <li>ii. Total no. of watch per box</li> </ol> </li> </ol> </li> </ol>			

TECHNICAL SPECIFICATION			
Republic of the Philippines Food and Drug Administration			
Lot No. 1 Item No.	Watch -Female for 35 years in government service	Qty. / Unit	7
Name of Manufacturer:		Country of Origin:	
Brand:		Model: (if applicable)	
Total ABC for Lot No. 1:			
PURCHASER'S SPECIFICATION		STATEMENT OF COMPLIANCE	
<p><b>A. Detailed Technical Specifications:</b></p> <ol style="list-style-type: none"> <li>1. Case size: at least 36mm</li> <li>2. Case thickness: 8-10mm</li> <li>3. Case material: Stainless steel</li> <li>4. Glass Material: Mineral Glass with sapphire</li> <li>5. Strap Materials: stainless steel; width: at least 16mm</li> <li>6. Dial, Case &amp; Strap Color: Must be complementary (Gold)</li> <li>7. Movement: Quartz Multifunction</li> <li>8. Water Resistant</li> <li>9. Scratch Resistant</li> <li>10. Watch Closure: Fold over clasp with push button</li> <li>11. Branded</li> <li>12. Warranty: 1 year</li> </ol>			
<p><b>B. Upon delivery the following shall be complied:</b></p> <ol style="list-style-type: none"> <li>1. Each watch must have an attached tag -- silver plated tag with embossed/engraved FDA logo and no. of years in service</li> <li>2.. Packaging Instructions: <ol style="list-style-type: none"> <li>a. Primary Packaging: Each piece is individually place in a color light gray PU leather box with pillow</li> <li>b. Secondary Packaging: individually packed, watched must be placed in a bigger box with the following label: <ol style="list-style-type: none"> <li>i. Watch for: specify gender and no. of years in service</li> <li>ii. Total no. of watch per box</li> </ol> </li> </ol> </li> </ol>			



**TECHNICAL SPECIFICATION**

Republic of the Philippines  
Food and Drug Administration

Lot No. 1	Watch -Male for 30 years in government service	Qty. / Unit	8
Item No.			

Name of Manufacturer:	Country of Origin:
-----------------------	--------------------

Brand:	Model: (if applicable)
--------	------------------------

Total ABC for Lot No. 1:

**PURCHASER'S SPECIFICATION****STATEMENT OF COMPLIANCE****A. Detailed Technical Specifications:**

1. Case size: at least 44mm
2. Case thickness: at least 8mm
3. Case material: Stainless steel
4. Glass Material: Mineral sapphire glass
5. Strap Materials: Silver stainless steel; width at least 22mm
6. Dial, Case & Strap Color: Must be complementary (Gold)
7. Movement: Quartz
8. Water Resistant
9. Watch Closure: 2 Pusher Fold over Clasp
10. Branded
11. Warranty: 1 year

**B. Upon delivery the following shall be complied:**

1. Each watch must have an attached tag -- silver plated tag with embossed/engraved FDA logo and no. of years in service
- 2.. Packaging Instructions:
  - a. Primary Packaging: Each piece is individually place in a color light gray PU leather box with pillow
  - b. Secondary Packaging: individually packed, watched must be placed in a bigger box with the following label:
    - i. Watch for: specify gender and no. of years in service
    - ii. Total no. of watch per box

TECHNICAL SPECIFICATION			
Republic of the Philippines Food and Drug Administration			
Lot No. 1 Item No.	Watch -Female for 30 years in government service	Qty. / Unit	21
Name of Manufacturer:		Country of Origin:	
Brand:		Model: (if applicable)	
Total ABC for Lot No. 1:			
PURCHASER'S SPECIFICATION		STATEMENT OF COMPLIANCE	
<b>A. Detailed Technical Specifications:</b>  1. Case size: at least 35mm 2. Case thickness: 8-10mm 3. Case material: Stainless steel 4. Glass Material: Mineral Glass with sapphire 5. Strap Materials: Silver stainless steel; width: at least 16mm 6. Dial, Case & Strap Color: Must be complementary (Gold) 7. Movement: Quartz Multifunction 8. Water Resistant 9. Scratch Resistant 10. Watch Closure: Double Press Deployant 11. Branded 12. Warranty: 1 year			
<b>B. Upon delivery the following shall be complied:</b>  1. Each watch must have an attached tag -- silver plated tag with embossed/engraved FDA logo and no. of years in service 2.. Packaging Instructions: a. Primary Packaging: Each piece is individually place in a color light gray PU leather box with pillow b. Secondary Packaging: individually packed, watched must be placed in a bigger box with the following label: i. Watch for: specify gender and no. of years in service ii. Total no. of watch per box			

TECHNICAL SPECIFICATION			
Republic of the Philippines Food and Drug Administration			
Lot No. 1 Item No.	Watch -Male for 25 years in government service	Qty. / Unit	3
Name of Manufacturer:		Country of Origin:	
Brand:		Model: (if applicable)	
Total ABC for Lot No. 1:			
PURCHASER'S SPECIFICATION		STATEMENT OF COMPLIANCE	
<b>A. Detailed Technical Specifications:</b>  1. Case size: at least 42mm 2. Case thickness: at least 11mm 3. Case material: Stainless steel 4. Glass Material: Mineral sapphire glass 5. Strap Materials: Silver stainless steel; width at least 22mm 6. Dial, Case & Strap Color: Must be complementary (white/silver/black) 7. Movement: Quartz 8. Water Resistant 9. Watch Closure: 2 Pusher Fold over Clasp 10. Branded 11. Warranty: 1 year			
<b>B. Upon delivery the following shall be complied:</b>  1. Each watch must have an attached tag -- silver plated tag with embossed/engraved FDA logo and no. of years in service 2.. Packaging Instructions: a. Primary Packaging: Each piece is individually place in a color light gray PU leather box with pillow b. Secondary Packaging: individually packed, watched must be placed in a bigger box with the following label: i. Watch for: specify gender and no. of years in service ii. Total no. of watch per box			

TECHNICAL SPECIFICATION			
Republic of the Philippines Food and Drug Administration			
Lot No. 1 Item No.	Watch -Female for 25 years in government service	Qty. / Unit	16
Name of Manufacturer:		Country of Origin:	
Brand:		Model: (if applicable)	
Total ABC for Lot No. 1:			
PURCHASER'S SPECIFICATION		STATEMENT OF COMPLIANCE	
<p><b>A. Detailed Technical Specifications:</b></p> <ol style="list-style-type: none"> <li>Case size: at least 36mm</li> <li>Case thickness: 6-8mm</li> <li>Case material: Stainless steel</li> <li>Glass Material: Mineral Glass</li> <li>Strap Materials: Silver stainless steel; width: at least 14mm</li> <li>Dial, Case &amp; Strap Color: Must be complementary (silver)</li> <li>Movement: Japanese Quartz</li> <li>Water Resistant</li> <li>Scratch Resistant</li> <li>Watch Closure: Double Press Deployant</li> <li>Branded</li> <li>Warranty: 1 year</li> </ol>			
<p><b>B. Upon delivery the following shall be complied:</b></p> <ol style="list-style-type: none"> <li>Each watch must have an attached tag -- silver plated tag with embossed/engraved FDA logo and no. of years in service</li> <li>Packaging Instructions: <ol style="list-style-type: none"> <li>Primary Packaging: Each piece is individually place in a color light gray PU leather box with pillow</li> <li>Secondary Packaging: individually packed, watched must be placed in a bigger box with the following label: <ol style="list-style-type: none"> <li>Watch for: specify gender and no. of years in service</li> <li>Total no. of watch per box</li> </ol> </li> </ol> </li> </ol>			



TECHNICAL SPECIFICATION			
Republic of the Philippines Food and Drug Administration			
Lot No. 1 Item No.	Watch -Male for 20 years in government service	Qty. / Unit	7
Name of Manufacturer:		Country of Origin:	
Brand:		Model: (if applicable)	
Total ABC for Lot No. 1:			
PURCHASER'S SPECIFICATION		STATEMENT OF COMPLIANCE	
<p><b>A. Detailed Technical Specifications:</b></p> <ol style="list-style-type: none"> <li>Case width: at least 40mm</li> <li>Case thickness: at least 10mm</li> <li>Case material: Stainless steel</li> <li>Glass Material: Mineral sapphire glass</li> <li>Strap Materials: Stainless steel; width at least 22mm</li> <li>Dial, Case &amp; Strap Color: Must be complementary (silver/black/white)</li> <li>Movement: Quartz</li> <li>Water Resistant</li> <li>Watch Closure: Single Pusher Fold over Clasp</li> <li>Branded</li> <li>Warranty: 1 year</li> </ol>			
<p><b>B. Upon delivery the following shall be complied:</b></p> <ol style="list-style-type: none"> <li>Each watch must have an attached tag -- silver plated tag with embossed/engraved FDA logo and no. of years in service</li> <li>Packaging Instructions: <ol style="list-style-type: none"> <li>Primary Packaging: Each piece is individually place in a color light gray PU leather box with pillow</li> <li>Secondary Packaging: individually packed, watched must be placed in a bigger box with the following label: <ol style="list-style-type: none"> <li>Watch for: specify gender and no. of years in service</li> <li>Total no. of watch per box</li> </ol> </li> </ol> </li> </ol>			

**TECHNICAL SPECIFICATION**

Republic of the Philippines  
Food and Drug Administration

Lot No. 1	Watch -Female for 20 years in government service	Qty. / Unit	12
Item No.			

Name of Manufacturer:

Country of Origin:

Brand:

Model: (if applicable)

Total ABC for Lot No. 1:

**PURCHASER'S SPECIFICATION****STATEMENT OF COMPLIANCE****A. Detailed Technical Specifications:**

1. Case size: at least 30mm
2. Case material: Stainless steel
3. Glass Material: Mineral Glass
4. Strap Materials: Silver stainless steel; alloy band
5. Dial, Case & Strap Color: Must be complementary (white/silver/black)
6. Movement: Japanese Quartz
7. Water Resistant
8. Scratch Resistant
9. Watch Closure: Clasp
10. Branded
11. Warranty: 1 year

**B. Upon delivery the following shall be complied:**

1. Each watch must have an attached tag -- silver plated tag with embossed/engraved FDA logo and no. of years in service
- 2.. Packaging Instructions:
  - a. Primary Packaging: Each piece is individually place in a color light gray PU leather box with pillow
  - b. Secondary Packaging: individually packed, watched must be placed in a bigger box with the following label:
    - i. Watch for: specify gender and no. of years in service
    - ii. Total no. of watch per box

TECHNICAL SPECIFICATION			
Republic of the Philippines Food and Drug Administration			
Lot No. 1 Item No.	Watch -Female for 15 years in government service	Qty. / Unit	6
Name of Manufacturer:		Country of Origin:	
Brand:		Model: (if applicable)	
Total ABC for Lot No. 1:			
PURCHASER'S SPECIFICATION		STATEMENT OF COMPLIANCE	
<b>A. Detailed Technical Specifications:</b>  1. Case size: at least 35mm 2. Case material: Stainless steel 3. Glass Material: Mineral Glass 4. Strap Materials: Leather; black 5. Dial, Case & Strap Color: Must be complementary (white/silver/black) 6. Movement: Quartz 7. Dial Hands: Basic 3 hands 8. Water Resistant 9. Scratch Resistant 10. Watch Closure: Buckle 11. Branded 12. Warranty: 1 year			
<b>B. Upon delivery the following shall be complied:</b>  1. Each watch must have an attached tag -- silver plated tag with embossed/engraved FDA logo and no. of years in service 2.. Packaging Instructions: a. Primary Packaging: Each piece is individually place in a color light gray PU leather box with pillow b. Secondary Packaging: individually packed, watched must be placed in a bigger box with the following label: i. Watch for: specify gender and no. of years in service ii. Total no. of watch per box			

**TECHNICAL SPECIFICATION**

Republic of the Philippines  
Food and Drug Administration

Lot No. 1	Watch -Male for 10 years in government service	Qty. / Unit	6
Item No.			

Name of Manufacturer:

Country of Origin:

Brand:

Model: (if applicable)

Total ABC for Lot No. 1:

**PURCHASER'S SPECIFICATION****STATEMENT OF COMPLIANCE****A. Detailed Technical Specifications:**

1. Case size: at least 36mm
2. Case material: Stainless steel
3. Glass Material: Mineral Glass
4. Strap Materials: Leather; black
5. Dial, Case & Strap Color: Must be complementary (white/silver/black)
6. Movement: Quartz
7. Water Resistant
8. Watch Closure: Pin buckle fastening
9. Branded
10. Warranty: 1 year

**B. Upon delivery the following shall be complied:**

1. Each watch must have an attached tag -- silver plated tag with embossed/engraved FDA logo and no. of years in service
- 2.. Packaging Instructions:
  - a. Primary Packaging: Each piece is individually place in a color light gray PU leather box with pillow
  - b. Secondary Packaging: individually packed, watched must be placed in a bigger box with the following label:
    - i. Watch for: specify gender and no. of years in service
    - ii. Total no. of watch per box



**TECHNICAL SPECIFICATION**

Republic of the Philippines  
Food and Drug Administration

Lot No. 1	Watch -Female for 10 years in government service	Qty. / Unit	12
Item No. 1			

**Name of Manufacturer:****Country of Origin:****Brand:****Model: (if applicable)****Total ABC for Lot No. 1:****PURCHASER'S SPECIFICATION****STATEMENT OF COMPLIANCE****A. Detailed Technical Specifications:**

1. Case size: at least 32mm
2. Case material: Stainless steel
3. Glass Material: Mineral Glass
4. Strap Materials: Leather; black
5. Dial, Case & Strap Color: Must be complementary (white/silver/black)
6. Movement: Quartz
7. Dial Hands: Basic 3 hands
8. Water Resistant
9. Scratch Resistant
10. Watch Closure: Buckle
11. Branded
12. Warranty: 1 year

**B. Upon delivery the following shall be complied:**

1. Each watch must have an attached tag -- silver plated tag with embossed/engraved FDA logo and no. of years in service
- 2.. Packaging Instructions:
  - a. Primary Packaging: Each piece is individually place in a color light gray PU leather box with pillow
  - b. Secondary Packaging: individually packed, watched must be placed in a bigger box with the following label:
    - i. Watch for: specify gender and no. of years in service
    - ii. Total no. of watch per box